### **READY TO USE**

PREVal

## Deliberately. Different

Prevail<sup>™</sup> Ready to Use Disinfectants clean by removing tough soils such as dirt, grime, mold, mildew, blood, urine, and fecal matter. Prevail<sup>™</sup> Ready to Use Disinfectants reduce odor-causing bacteria leaving surfaces smelling clean and fresh.

Its non-abrasive formula is designed for use on hard, nonporous surfaces associated with non-food contact countertops, cages, kennels and animal equipment and devices found in animal housing facilities.

#### **AVAILABLE PRODUCTS**

Prevail<sup>™</sup> Animal Premise Disinfectant Cleaner & Deodorizer Ready to Use

#### DIN: 02436795

SHELF LIFE: 2 Years

**CLAIMS:** 3 Minutes (Tuberculocidal, Virucidal, Bactericidal, Fungicidal) 30 Seconds (Broad-Spectrum Sanitizer)

1. Format Size: 1 L Bottle Product Code: 12101 Packaging Specifications: 12x1 L

2. Format Size: 4 L Bottle Product Code: 12105 Packaging Specifications: 4x4 L

For more information, visit: ViroxAnimalHealth.com

# Animal Premise Disinfectant Cleaner & Deodorizer Ready to Use formulated with

**Faster** Effective against a broad spectrum of pathogens in only 3 minutes, including Canine parvovirus

**Cleaner** One-Step Disinfectant Cleaners with powerful surfactants (detergents) simplifying sanitation of animal environments saving time and money

**Responsible** Prevail<sup>™</sup> is essentially non-toxic<sup>\*</sup> and non-irritating, does not require the use of PPE

Sustainable Hydrogen peroxide breaks down into water and oxygen

**Compatible** Designed to be compatible with a wide range of materials, protecting your assets

Prevail<sup>™</sup> was designed to be easier on both animals and humans improving user compliance. The ready to use formats require no mixing or measuring, improving protocol compliance.

Prevail<sup>™</sup> Disinfectants are ideal for cleaning and disinfection of clinical areas and surfaces such as work stations, fast disinfection of non-critical devices and equipment, and maintaining a clean environment to reduce the potential for cross-contamination or the spread of infections.

#### CLEANING, DISINFECTION AND DEODORIZING OF EQUIPMENT AND NON-POROUS HARD SURFACES

**1.** Pre-clean heavily soiled areas.

- 2. Apply directly onto surface until thoroughly wet.
- 3. Ensure surface remains wet for 3 minutes.

\* Lowest WHMIS toxicity category for all 3 criteria: oral, dermal and inhalation \*\* Refer to reference sheet for additional claims and use instructions



2770 Coventry Rd., Oakville, ON L6H 6R1 1-800-387-7578 | viroxanimalhealth.com





#### FOR USE IN...

- Agricultural farms and barns
- Animal housing facilities and areas
- Animal shelters and kennels
- Barn Premise
- Farm Premise
- Pet housing facilities
- Surgical suites
- Transportation vehicles
- Treatment facilities
- Veterinary practices and offices

#### FOR USE ON...

#### Hard Surfaces such as:

- Counters
- Desks
- Cage spot cleaning
- Door & drawer handles
- Light switches
- Faucets

#### Devices & Equipment such as:

- Stethoscopes
- Otoscope handles
- Clipper handles
- Microscope knobs
- Tack
- Hose ends
- Computer keyboards & mice
- Phones

#### **CERTIFICATIONS**

EcoLogo/UL Approved