

FOR VETERINARY USE

PREvail®



**Deliberately
Different.**



For use in:

- Animal housing facilities, shelters & kennels
- Grooming facilities
- Pet stores
- Surgical suites & treatment facilities
- Transportation vehicles
- Veterinary practices & offices
- Other small animal facilities



For use on:

- Animal crates & carts
- Automatic feeders
- Counters, exam tables & other hard surfaces
- Diagnostic & laboratory equipment
- Dog runs
- Exam room instruments
- Food & water bowls
- Grooming & handling equipment
- Harnesses & leashes
- Transportation crates & vehicles

Animal Premise Disinfectant Cleaner & Deodorizer Concentrate formulated with



Prevail™ Disinfectants are Deliberately Different™ from other legacy disinfectant products available today. Our unique and innovative Ready to Dispense (RTD¹) format was designed to make disinfection *easy, reliable, accurate* and *gentle* for people and animals.

Cleaner Prevail™ Concentrate has superior cleaning efficacy with no harsh chemicals or offensive odours. Prevail™ Concentrate is an effective cleaner and deodorizer eliminating the smell of feces, pheromones, urine and other soil odours. The smell of clean is no smell.

Faster Prevail™ Concentrate provides confidence that disinfection will occur. Prevail™ Concentrate was designed to be effective against a broad spectrum of pathogens including Canine parvovirus in only 5 minutes, improving turnover times and reducing the risk of outbreaks.

Responsible Prevail™ Concentrate RTD does not require the use of PPE and is gentle for staff, animals, and the environment. Prevail™ Concentrate was designed without the use of fragrances and does not have strong chemicals odours, improving indoor air quality and reducing fear, anxiety and stress on animals.

Sustainable Prevail™ Concentrate provides assurance that the environment will not be negatively impacted. The active ingredient, hydrogen peroxide, breaks down into water and oxygen, leaving no active residues behind and does not require any special disposal methods.

Compatible Prevail™ Concentrate is compatible with a wide range of materials. Prevail™ Concentrate will protect your assets by reducing corrosion and wear.

* See reference sheet for a complete list of pathogenic organisms, additional features, claims, and directions for use.

Why choose Ready to Dispense (RTD[†])?

- Ready to Dispense offers a Deliberately Different™ method to simply, accurately, reliably, and safely dispense a concentrate into a bottle or bucket without the need for dilution equipment.
- **Comfort & Confidence:**
 - The ready to dispense technology ensures the right concentration of Prevail™ every time, reducing misuse, spillage or overuse resulting in cost and time savings.
 - The easy to use head was designed for both bottle and bucket fill by connecting the bottle to a water source, ideal for areas with poor or inconsistent water pressure.
 - The built in dilution head eliminates exposure to the concentrate solution, reducing the need for PPE and minimizing the need for staff and volunteer training.

Four Easy Steps to Disinfection

1 Open  Open the case and read the instructions for use.

2 Connect  Connect Prevail™ RTD to a water supply.

3 Select  Select the desired dilution and fill rate with easy to use icons.

4 Fill  Squeeze handle and fill.
You're ready to disinfect!

For more information, visit:
ViroxAnimalHealth.com

2770 Coventry Rd., Oakville, ON L6H 6R1 CANADA
1-800-387-7578 | viroxanimalhealth.com



Prevail™ helps to ensure a Fear Free™ environment by eliminating harsh chemical odours that trigger fear, anxiety and stress for pets and pet owners during veterinary visits. Virox Animal Health™ is a proud partner of Fear Free™, developed by Dr. Marty Becker, VetFolio and the American Animal Hospital Association.

Available Products

Prevail™ Animal Premise Disinfectant Cleaner & Deodorizer

DIN 02436809

- 1. Format Size:** 5L RTD
- Product Code:** 12308
- Packaging Specifications:** 5L RTD

Other Formats Available

Claims

5 Minutes @ 1:40 Dilution
(*Virucidal, Bactericidal, Fungicidal, Tuberculocidal)

30 Seconds @ 1:128 Dilution
(Broad-Spectrum Sanitizer)

Shelf Life

2 Years, 30 Days Diluted (meets stability testing for up to 30 days after dilution with tap water)

[†]Dilution control and dispensing technology utilized under license from Diversy Care, a division of Sealed Air Corporation.